# MINILITH® SL1

SPEICHELSTEINE NICHTINVASIV BEHANDELN  
THE NON-INVASIVE THERAPIE OF SALIVARY STONES

STORZ MEDICAL

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## MINILITH® SL1 THE INNOVATIVE SOLUTION

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The MINILITH® SL1 from STORZ MEDICAL is the first lithotripter worldwide designed for the extracorporeal treatment of salivary stones.

The MINILITH® SL1 is specially adapted for the gentle and efficient treatment of stones in the submandibular, parotid and sublingual glands and is a proven alternative to surgery. The treatment

- does not require an anesthetic
- can be performed on an outpatient basis and
- is gentle and efficient.

The MINILITH® SL1 consists of an electromagnetic pressure wave source, the core of the lithotripter, an in-line ultrasound system and a supply and control unit on wheels.

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Extracorporeal lithotripsy by means of shock or pressure waves was first carried out in human beings in 1980 and, within a few years, has revolutionized the treatment of urological stones. Five more years of development were needed before the technique could be used to treat gallstones. The lithotripters now available differ with regard to concretment localization, positioning of the patient and monitoring of stone disintegration.

Depending on the application of the lithotripter, the most important factors are the selection of an appropriate system for stone localization and of an appropriate energy source, the core component of any lithotripter.

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### **EXTRACORPOREAL LITHOTRIPSY OF SALIVARY STONES**

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For some years now it has been possible to use extracorporeal outpatient lithotripsy to fragment salivary stones to such an extent that they can be washed out with the saliva as sand-like material. Thanks to this method, surgical removal of the salivary gland can be avoided in most cases. In the past, surgical treatment was necessary when dilation or incision of the salivary duct failed to allow removal of the stone from within the oral cavity. The use of extracorporeal lithotripsy avoids not only the risks associated with a general anesthesia but also those of damage to the facial nerve and postoperative Frey syndrome.

Non-invasive stone disintegration is possible thanks to the precise localization provided by the coaxial integration of the in-line ultrasound transducer in the therapy source of the MINILITH® SL1. Because of the limited volume of the therapy focus, the surrounding tissue, the jawbones and the teeth are almost completely unaffected by the pressure waves.

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## LOCALIZATION OF SALIVARY STONES WITH THE MINILITH® SL 1

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The ultrasound localization system of the MINILITH® SL1 is *in-line* integrated. The scan plane of the ultrasound transducer can be rotated. This ensures optimal and precise localization and fragmentation control. The ultrasound transducer can be removed without difficulty from the treatment source, thereby allowing exact and comfortable pre-localization of the stone as well as general ultrasound diagnosis.

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
The treatment source of the MINILITH® SL1 is designed around the proven STORZ pressure wave principle and consists of an electromagnetic cylindrical coil with a parabolic reflector as the focusing unit. The focus depth, focus dimensions, size of the applicator and energy levels of the treatment source have been adjusted to the requirements of salivary stone lithotripsy, thus eliminating the need for a local anesthetic. After fragmentation, sand-like material is passed with the medically induced increase in saliva flow.

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
The treatment source is attached to a mobile arm and can be moved by hand. The arm has five axes of movement, each of which can be hydraulically locked and fixed. Fine regulation during the locked condition allows optimal adjustment to the shape of the patient's body. The pressure wave can be triggered by a hand or foot-operated switch as required. Existing couches, dentistry chairs or similar may be used for treatment during which the patient lies comfortably on his or her back.

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Clear image of stone near the submandibular gland.

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After treatment the patient is »stone-free«.

## STORZ MEDICAL - IHR SPEZIALIST IN LITHOTRIPSIE

## STORZ MEDICAL - YOUR SPECIALIST IN LITHOTRIPSY

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STORZ lithotripters have been successfully used all over the world for years. The principle of the cylindrical electromagnetic energy source has proved itself in practice again and again. Performance, efficiency and reliability are just some of the features that make the STORZ lithotripters unique in their field.

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With the MINILITH® SL1, STORZ MEDICAL has adjusted its tried and tested technology to meet the demands of salivary stone treatment. The results show that the MINILITH® SL1 is an important milestone in the gentle extracorporeal treatment of salivary stones and that it eases the burden on both patient and doctor.

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
The compact ergonomic system requires little space. Just plug the device in and it is ready for use. High-quality components guarantee a high degree of reliability and a long operating life.

Our worldwide service network guarantees comprehensive advice and customer service. Installation of the equipment and regular servicing are an integral part of our package.

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## TECHNICAL DATA

CE in accordance with the MDD directive (IEC 601-1 et seq)

<b>Dimensions</b>	(w x l x h)
Transport position 1	650 mm x 990 mm x 1480 mm
Transport position 2	706 mm x 1118 mm x 1307 mm
Weight (without ultrasound unit)	approx. 190 kg
Operating voltage	100 / 115 / 230 VAC, 50 / 60 Hz
Power consumption	max. 2.0 kVA (incl. external devices)
Closed water circuit with integrated water conditioning	

Therapy arm with three joints; freely moveable in three directions	
max. diameter of therapy head	155 mm
Range of action	1168 mm
Angular range	max. 330°
Height adjustment	550 mm - 1250 mm
Precise adjustment along two orthogonal axes (optional)	

Adjustable focus pressure* P <sub>r</sub>	12.5 - 91.8 MPa
Energy flux density* ED	0.15 - 1.49 mJ/mm <sup>2</sup>
Frequency	1 - 4 Hz
Focal size* f <sub>x</sub> (-6dB) x f <sub>z</sub> (-6 dB) (Level 6)	3.4 mm x 25 mm
Focal distance	0 - 50 mm
Introduction of energy over a large surface area	apertur angle 84°

### Ultrasound localization

7.5 MHz sector scanner, axially integrated into the shock wave source, the scanning level can be pivoted by 360°. The ultrasound unit can be used independently of the MINILITH® SL1.

### X-ray localization

Optional in-line targeting device.

\* measured with fiber optic hydrophone

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## **THERE IS MORE TO CUSTOMER SATISFACTION THAN JUST A GOOD PRODUCT**

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Even the best system and the best technique require regular maintenance. Our team of application engineers provide training for your personnel on your premises in order to make them familiar with this new type of therapy. It goes without saying that we will also be at your disposal at any later date, if need be. The components of the MINILITH® SL1 are easy to service. However, even if you require specialized help, you can rely on your local STORZ service centre to rectify errors within the shortest of time.

**STORZ MEDICAL AG · UNTERSEESTRASSE 47 · CH-8280 KREUZLINGEN  
TELEFON +41 (0)71 677 45 45 · TELEFAX +41 (0)71 677 45 05  
E-MAIL: [INFO@STORZMEDICAL.COM](mailto:INFO@STORZMEDICAL.COM) · [WWW.STORZMEDICAL.COM](http://WWW.STORZMEDICAL.COM)**

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